

# Cavity ringdown spectroscopy

What if we do absorption spectroscopy over a very  
long range?

What if the absorption feature we need to measure is really weak or if you need very high precision?

- What do you think we can do?

What if the absorption feature we need to measure is really weak, or if you need very high precision?

- Integrate longer
- Use a more sensitive detector
- Cool the detector to lower thermal noise in the system
- Control laser frequency to high degree
- **USE A REALLY LONG PATH SO THAT THE ABSORPTION IS GREATER**
  - But what's the limit of that? Can I make a path that's a kilometer long?


# Cavities! (not in your teeth hopefully)

- Two mirrors properly aligned such that they light is an incredibly high number of “bounces” in the mirror.
  - Gives a very long effective path length
  - High intensity in the laser cavity

# But how does it work?


- A laser pulse bounces around the cavity and a little energy from it leaks out each round trip
- You can approximate the peaks of each bounce together as an exponential decay
- Light intensity is an exponential function of extinction coefficient in the cavity

$$\begin{aligned}
 I_{\text{out}} &= I_{\text{in}} \exp \left[ - \left\{ \frac{(1 - R)c}{L} t + \left( \text{absorption loss for} \right) \left( \frac{\text{total number of}}{\text{sample per round trip}} \right) \left( \frac{\text{round trips}}{\text{round trips}} \right) \right\} \right] \\
 &= I_{\text{in}} \exp \left[ - \left\{ \frac{(1 - R)c}{L} t + (2\alpha d) \left( \frac{tc}{2L} \right) \right\} \right] \\
 &= I_{\text{in}} \exp \left[ - \frac{tc}{L} \{ (1 - R) + (\alpha d) \} \right] \\
 &= I_{\text{in}} \exp \left[ - \frac{t}{\tau} \right]; \quad \text{where } \tau = \frac{L}{[(1 - R) + \alpha d]c}
 \end{aligned}$$



# How do we actually get the data from that exponential0

- Take the natural log of the exponential decay
- Fit a line to it
- Solve for the absorption



# Problems

- Very narrowband measurements
- Slow since you need to record multiple ringdowns
- **VERY FINNICKY TO SETUP**
  - You have at least four mirrors that need to be set up just-so to actually get it where you need it
- Requires some way to modulate the cavity or power within it
  - Pulsed laser
    - Diode lasers are great for this because optical feedback into it, can lock it to a frequency; See “external cavity diode lasers”
  - Modulate cavity off resonance
- Need a fast detector
- Need some really good, high reflectance mirrors
  - Normally we have 99.9% reflective , to do this well you need 99.999% reflective
- Need to match the lasers’ spatial modes to the spatial modes of the cavity

# What You're gonna build

